


Calibration Slides on Glass slides 76 x25mm so you can set your spacial calibration in software or eyepiece graticules.


TS-TCM-H With scales 10mm in 100 div. , 1mm in 100 div , dot 0.15mm , dot 0.07mm



TS-TCM-L Standard slide 1mm in 100 div each div 0.01mm



TS-TCM-S with scale 10mm in 100 div ideal for stereo and macro microscopy



Supplied with plastic case.



Scientific Instrument and Optical Sales
Unit 1, 62 Bishop Street
Kelvin Grove, Brisbane,
QLD Australia 4059

Ph : (07) 3356 0233
Fax : (07) 3356 0344
email: sios@pacific.net.au